


Exploring the iGEM Distribution



Overview

Getting Started

Introduction, information and structure

Data Exploration & Analysis

Diving deep into the Registry, analyzing how the distribution is used by teams

Suggestions

Files and software demo

Discussion

Registry of Standard Biological Parts

[tools](#) [catalog](#) [repository](#) [assembly](#) [protocols](#) [help](#) [search](#)

BBa_



DNA Repository Plates and Boxes

Physical DNA is held in tubes in freezer boxes or multi-well plates. This program manages the contents of boxes and plates.

Label: **2021 Kit Plate 1** ID: 6411
Description: Spring 2021 DNA Distribution Kit Plate 1 384-Well Plate
Location: 2021-07-02 10:59:48
Substance: DNA
Plate Status: OK
Aliquot: ☐ Checked if this plate contains unprocessed samples from its source plate.

[Get antibiotic files for this plate](#)
[Get an Excel file for this plate](#)
[Get a detailed Excel file for this plate](#)

Gel Images and Results

[Wells 1A thru 6H](#)
[Wells 7A thru 12H](#)

Plate Images and Results

[Add a plate image](#)

Sequencing and Results

[Go to sequencing](#)

Contents:

Well	Part	Plasmid	Resist. Cell	Comments
Quality control information: Sequencing, Antibiotics, Restriction Digests				
1A	BBa_K314110	pSB1C3	E. coli strain NEB 10-beta	
1B	BBa_K731722	pSB1C3	E. coli strain NEB 10-beta	
1C	BBa_K398326	pSB1C3	E. coli strain NEB 10-beta#1 - 520 ng	
1D	BBa_K731722	pSB1C3	E. coli strain NEB 10-beta	I was not able to remove this line from the list. It is the same sample as #21
1E	BBa_K398331	pSB1C3	E. coli strain NEB 10-beta#2 - 990 ng	
1F	BBa_K808025	pSB1C3	E. coli strain NEB 10-beta	may be toxic for gramm negativ bacteria due to lipase activity if its expressed into periplasma or surface expressed
1G	BBa_K314100	pSB1C3	E. coli strain NEB 10-beta	
1H	BBa_K808001	pSB1C3	E. coli strain NEB 10-beta	
1I	BBa_K314101	pSB1C3	E. coli strain NEB 10-beta	
1J	BBa_K808003	pSB1C3	E. coli strain NEB 10-beta	
1K	BBa_K314201	pSB1C3	E. coli strain NEB 10-beta	
1L	BBa_K808010	pSB1C3	E. coli strain NEB 10-beta	
1M	BBa_K314202	pSB1C3	E. coli strain NEB 10-beta	
1N	BBa_K808011	pSB1C3	E. coli strain NEB 10-beta	
1O	BBa_K346002	pSB1C3	E. coli strain NEB 10-beta	PmerT promoter
1P	BBa_K808013	pSB1C3	E. coli strain NEB 10-beta	
2A	BBa_K548000	pSB1C3	E. coli strain NEB 10-beta	
2B	BBa_K325210	pSB1C3	E. coli strain NEB 10-beta	

Parsing the data

Files are not Excel but CSV files

Some software applications like Excel and Pandas have problems parsing the file even though the file seems to use RFC4180

Deleting the space between commas and double quotes can easily fix this problem

```
----> 1 test_df = pd.read_csv('old_formatting.csv')

c:\users\cedkb\appdata\local\programs\python\python37\lib
p, delimiter, header, names, index_col, usecols, squeeze,
se_values, skipinitialspace, skiprows, skipfooter, nrows,
se_dates, infer_datetime_format, keep_date_col, date_pars
s, decimal, lineterminator, quotechar, quoting, doublequ
lines, delim_whitespace, low_memory, memory_map, float_pr
686     )
687
--> 688     return _read(filepath_or_buffer, kwds)
689
690

c:\users\cedkb\appdata\local\programs\python\python37\lib
458
459     try:
--> 460         data = parser.read(nrows)
461     finally:
462         parser.close()

c:\users\cedkb\appdata\local\programs\python\python37\lib
1196     def read(self, nrows=None):
1197         nrows = _validate_integer("nrows", nrows,
-> 1198         ret = self._engine.read(nrows)
1199
1200         # May alter columns / col_dict

c:\users\cedkb\appdata\local\programs\python\python37\lib
2155     def read(self, nrows=None):
2156         try:
-> 2157         data = self._reader.read(nrows)
2158         except StopIteration:
2159             if self._first_chunk:

pandas\_libs\parsers.pyx in pandas._libs.parsers.TextReader
pandas\_libs\parsers.pyx in pandas._libs.parsers.TextReader
pandas\_libs\parsers.pyx in pandas._libs.parsers.TextReader
pandas\_libs\parsers.pyx in pandas._libs.parsers.TextReader
pandas\_libs\parsers.pyx in pandas._libs.parsers.raise_p
ParserError: Error tokenizing data. C error: Expected 17
```

Do the csv files match the data on the website?

All ccdB parts have been withdrawn from the Registry. Samples of parts containing the ccdB gene cannot be requested. - iGEM HQ

They do in fact differ slightly. The csv file contains more entries than the website

The parts in question are:

BBa_K581008

BBa_K145151

BBa_K805013

BBa_K658001

All of these include the ccdB gene which was removed from the registry due to patent issues

CSV overview

Duplicate Columns:

Plasmid

Empty Columns:

Resistance, Gel Overall, Quantity, Seq Comment








Just one unique value:


Sequencing, Well Status

Columns

Plate
Well
Part
Plasmid
Well Status
Comments
Type
Subparts
Source
Gel Overall
Quantity
Plasmid
Sequencing
Seq Comment
Short Description


Sequencing

-  **Confirmed**
-  **Partially Confirmed** - software is only able to partially confirm sequence, most likely due to one read being poor
-  **Long Part** – length of sequence reads are insufficient to cover the middle of the part
-  **Inconsistent** - part does not match its target sequence, may have a single bp mutation or not match at all
-  **Bad Sequence** – usually caused by low DNA concentration or incorrect primers
-  **Single Error**
-  **No Sequencing Information**




Sequencing II

Differentiate between
confirmed and all others
categories




Type of parts

Coding	577
Composite	329
Regulatory	254
Reporter	239
Generator	126
Intermediate	106
Translational_Unit	82
Measurement	76
DNA	66
Device	63
Terminator	32
Signalling	32
RBS	27
Protein_Domain	26
RNA	24
Project	19
Other	15
Tag	13
Inverter	8
Conjugation	4
T7	3
Temporary	1



Plasmid


pSB1C3	1936
pSB1A2	83
pSB1C6	34
BBa_J61002	16
pSB1C00	9
pSB3K3	5
pSB1A3	3
pSB1AK3	2
BBa_J63009	2
pSB1AK8	2
BBa_P10501	1
pSB2K3	1
pSB1K01	1
BBa_P10506	1
pSB4K5	1
pSB1K3	1
pSB3C04	1
pSB1K04	1
pSB6A1	1
pSB3C02	1
pSB1AT3	1
pSB1K03	1
pSB3C01	1
BBa_J63010	1
BBa_P10503	1
BBa_P10504	1
BBa_P10507	1
pSB1K02	1
pSB1T3	1
pSB1AC3	1
BBa_P10502	1
pSB4C5	1
pSB3C5	1
pSB1A10	1
pSB3T5	1
BBa_P10505	1
pSB4A5	1
BBa_J23006	1
BBa_P10508	1
pSB3C03	1



GoldenBraid destination vector pDGB3 alpha1R

More information to come

Sequence and Features



Note: The sequence below is incorrect because the backbone that the SP device was inserted into was psb1A2 rather than psb1A3, so there is a terminator (not annotated) in the sequence below at position 1044-1121 that isn't actually present in the DNA of psb1A10. Updating this will take some time because I need to change all the feature locations by hand - but the [Vector NTI file](#) is correct.

Some of the plasmids that are used lack basic data


Some of the plasmids used in the distribution plate lack either:

- General information
- Information about the resistance
- Sequence annotations

Or everything of the above.

Additionally, some of the information may also be incorrect

Registry of Standard Biological Parts


[tools](#)
[catalog](#)
[repository](#)
[assembly](#)
[protocols](#)
[help](#)
[search](#)

BBa_



Repository: Eppendorf Files for Broth Plates

[Show all plates](#)

Label:	2021 Kit Plate 1	ID: 6411
Description:	Spring 2021 DNA Distribution Kit Plate 1	384-Well Plate
Location:		2021-07-02 10:59:48
Substance:	DNA	

Well	Part	Plasmid	Cell
Unexpected Plasmid: pSB1C3 in this well:			
1A	BBa_K314110	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1B	BBa_K731722	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1C	BBa_K398326	pSB1C3	E. coli strain NEB 10-beta #1 - 520 ng
Unexpected Plasmid: pSB1C3 in this well:			
1D	BBa_K731722	pSB1C3	E. coli strain NEB 10-beta I was not able to remove this line from the list. It is the same sample as #21
Unexpected Plasmid: pSB1C3 in this well:			
1E	BBa_K398331	pSB1C3	E. coli strain NEB 10-beta #2 - 990 ng
Unexpected Plasmid: pSB1C3 in this well:			
1F	BBa_K808025	pSB1C3	E. coli strain NEB 10-beta may be toxic for gramm negativ bacteria due to lipase activity if its expressed into periplasma or surface expressed
Unexpected Plasmid: pSB1C3 in this well:			
1G	BBa_K314100	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1H	BBa_K808001	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1I	BBa_K314101	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1J	BBa_K808003	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1K	BBa_K314201	pSB1C3	E. coli strain NEB 10-beta
Unexpected Plasmid: pSB1C3 in this well:			
1L	BBa_K808016	pSB1C3	E. coli strain NEB 10-beta

Schrödinger Parts

Some of the comments suggest that the part is in a different plasmid than officially stated.

Comments tend to add uncertainty rather than provide helpful insights.



10H	BBa_E0020	pSB1C00	pSB1C5 is a Type IIS Plasmid for Level 0 parts
10J	BBa_E0030	pSB1C00	pSB1C5 is a Type IIS Plasmid for Level 0 parts
10L	BBa_J97004	pSB1C00	pSB1C5 is a Type IIS Plasmid for Level 0 parts
10N	BBa_J97005	pSB1C00	pSB1C5 is a Type IIS Plasmid for Level 0 parts
10P	BBa_J97006	pSB1C00	pSB1C5 is a Type IIS Plasmid for Level 0 parts

15D	BBa_I14044	pSB1C3	E. coli strain NEB 10-betaPlasmid unknown, AmpR
15E	BBa_K863020	pSB1C3	E. coli strain NEB 10-beta
15F	BBa_I13013	pSB1C3	E. coli strain NEB 10-betaPlasmid unknown, AmpR

13P [BBa_K1093002](#) pSB1C3

Sorry for the Plasmid Backbone being pSB1A2 because of having no time to transfer.

404

File not found

What about files?

A *in silico* version of the distribution (i.e. GenBank files) would be helpful for both exploring as well as cloning parts of interest

Sadly, iGEM does not offer the possibility to download such files

Furthermore, the linked part tool with which one could download a part as a genbank file is no longer accessible

Plasmid Mapper


Although the registry links to a plasmid mapper, it is not a substitute for raw data

It doesn't feature the actual plasmid maps, it only creates a custom in silico cloning of the part

Contains less than half the plasmids used in the distribution

It's currently broken for all major browsers





Relationship between parts in the distribution

Embedded all of the parts into a
directed graph


The resulting graph consists of 1,985
edges and 10,236 vertices

Graph

To better visualize the graph we utilize the force-directed graph algorithm from Fruchterman & Reingold

The resulting graph is surprisingly connected

One part in particular stands out, it seems to be included in every other part



Mysterious Part

A rather specialized part without defined uses

Only consist of 1bp... A

This false sequence occurrence both on the registry site as well as the registry database dump

Part:BBa_K1497024

Designed by: Sascha Hein, Rene Sahm Group: iGEM14_TU_Darmstadt (2014-10-07)

GBD-Domain

The GBD domain and its ligand (BBa_K771106) are suitable tools for protein colocalization. Initially, the domain was a part of the N-WASP protein (GTPase binding domain) in *Rattus rattus*. The GBD domain is used as a binding unit of the so-called protein scaffold published by Dueber et al. in 2009. The scaffold (BBa_K1497033) is composed of different binding units, which enable the assembly of multiple target proteins. The Domain was initially edited by Team BioX-Shanghai 2012 (BBa_K771105). The iGEM Team TU Darmstadt 2014 modified this BioBrick by adding a BglII and BamHI restriction site in front of and behind the previously constructed domain sequence and codon optimized it for expression in *E. coli*.

Now, different binding units of the scaffold protein can be fused together without the introduction of restriction sites. This allows the easy construction of BioBricks of different permutations of the scaffold protein domains.





Figure 1 Crystal structure of the GBD domain from the protein (red) of *Rattus rattus* binding the small molecule (1-(3,6-Dibromo-carbazol-9-yl)-3-dimethylamino-propan-amine). The GBD domain is locked by Wiskostatin in its autoinhibited conformation (Peterson et al. 2004). PDB entry 1T84.

Subparts Ruler SS DS				Scars: Show Hide		Vertical: Show Hide		Length: 1 bp	
1	1 a	11	21	31	41	51	61	71	81

Assembly Compatibility: 10 12 21 23 25 1000

Duplicate Wells

<u>Part</u>	<u>Count</u>	<u>Name</u>	<u>Uses</u>
BBa_J04450	40	RFP Coding Device	69
BBa_E0240	5	GFP Generator	293
BBa_R0040	4	TetR Repressible Promoter	1081
BBa_P10599	4	GoldenBraid destination vector alpha insert	0
BBa_P10598	4	GoldenBraid destination vector omega insert	0
BBa_J04455	4	RFP Selection Device for Even Level Loop Type IIS Parts	0
BBa_J04454	4	RFP Selection Device for Odd Level Loop Type IIS Parts	0
BBa_E0020	3	Cyan Fluorescent Protein	95
BBa_I20270	3	Promoter MeasKit	2
BBa_J364000	3	Test Device 1 for the iGEM InterLab Study	8



Parts w/ same sequence but different name

amilCP blue/purple chromoprotein (incl RBS)

['BBa_K592025', 'BBa_K1033930']

Screening plasmid intermediate

['BBa_K1357010', 'BBa_I13507']

Promoter (luxR & HSL regulated -- lux pL)

['BBa_R0063', 'BBa_K783024']

green fluorescent protein derived from jellyfish

['BBa_K895006', 'BBa_E0040']

CFP

['BBa_E0020', 'BBa_K1418040']

RBS30

['BBa_B0030', 'BBa_K1789002']

RFP Selection Device for Level 0 Type IIS Parts

['BBa_J04450', 'BBa_J04452', 'BBa_J04455',
'BBa_J04454']

dCBD with N-terminal linker in RFC 25

['BBa_K1321340', 'BBa_K1321347']

MoClo format of BBa_J23112 with AB fusion

['BBa_K1114003', 'BBa_K1114014']

pTetR-lasR-Term-Term

['BBa_K876015', 'BBa_K876057']

pSB-AraC-pBAD

['BBa_I0500', 'BBa_K1321333']

constitutive promoter family member

['BBa_K823007', 'BBa_J23103', 'BBa_J23112']

['BBa_K823010', 'BBa_J23113']

['BBa_K823013', 'BBa_J23117']

['BBa_K823006', 'BBa_J23102']

['BBa_K823004', 'BBa_J23100']

['BBa_K823014', 'BBa_J23118']

['BBa_K823005', 'BBa_J23101']

['BBa_K783028', 'BBa_J23110']

['BBa_K823008', 'BBa_J23106']

['BBa_K783032', 'BBa_J23115']

['BBa_K823011', 'BBa_K783031', 'BBa_J23114']

Twin Parts

"Two or more parts are twins if they have the same DNA sequence"

Although twin parts should generally be avoided, there are parts that have over 50 twins

This makes it more difficult to collect all of the necessary data

```
▼ <twins>
  <twin>BBa_J34805</twin>
  <twin>BBa_J70612</twin>
  <twin>BBa_J72049</twin>
  <twin>BBa_K783020</twin>
  <twin>BBa_M36912</twin>
  <twin>BBa_K1045009</twin>
  <twin>BBa_T6000</twin>
  <twin>BBa_K1520029</twin>
  <twin>BBa_K1685000</twin>
  <twin>BBa_K2042020</twin>
  <twin>BBa_K1972011</twin>
  <twin>BBa_K2023000</twin>
  <twin>BBa_S05347</twin>
  <twin>BBa_K1903000</twin>
  <twin>BBa_K2175002</twin>
  <twin>BBa_K2101000</twin>
  <twin>BBa_K2406215</twin>
  <twin>BBa_K2560035</twin>
  <twin>BBa_K3136005</twin>
  <twin>BBa_K3228005</twin>
  <twin>BBa_K3296007</twin>
  <twin>BBa_K3158004</twin>
  <twin>BBa_K3135007</twin>
  <twin>BBa_K2969993</twin>
</twins>
```

T7 Promoter
 T7 promoter (strong promoter from T7 bacteriophage)
 UNS 3 Sequence, from Torella et al., 2013
 hPCD
 flyPCD
 fshPCD
 UNS 2 Sequence, from Torella et al., 2013
 yeast GAL1 promoter
 Medium strength T7.2 RBS
 T7 consensus -10 and rest
 QPI (B0034.C0051.B0015.R0051)
 QPI (B0034.C0040.B0015.R0040)
 TetR repressed POPS/RIPS generator
 WM_Pad1
 WM_Pad2
 Strong T7.2 RBS
 ADH1 terminator from S. cerevisiae
 Linker b (BamH I)
 QPI (B0034.C0012.B0015.R0011)
 SH3 domain + Linker
 PDZ domain + Linker
 Linker h
 ECK120029600 - Escherichia coli K-12 terminator
 His tag
 attB1 recombination site
 sfGFP
 VP64
 Terminator (artificial, small, %T~85%)
 MoClo RBS B0034
 [rnpB-T1] Terminator
 N-terminal start overhang (T)(A)-(G)ATG=RBS+Start RFC[105] A
 GST
 Linker e（cutting site of Prescission Protease）
 hixC binding site for Salmonella typhimurium Hin recombinase
 T7 R0.3 RNaseIII site
 Transcription Terminator (Strong)
 scaffold of sgRNA in CRISPR/Cas system
 Linker a（Nde I and Nhe I）
 Sumo tag
 Terminator (His)
 Stop codon free GFP in RFC[23] standard
 hU6 Promoter

Most used parts that are not in the distribution

Analyzed the usage of parts in both the distribution and the entire registry

While some of the most commonly used parts are missing, others are included in the registry under a different name which can further complicate their use



Most parts are barely used

Including all of the twin parts, the majority of parts in the distribution are rarely used


818 (41.8%) haven't been used at all

1633 (83.5%) have been used less than 10 times


Total Uses




Age of the Parts in the distribution




Number of Parts with over 50 uses




Adding up the uses for each year



Number of Parts **vs.** the amount of Subparts they consists of



Number of Parts with over 10 uses **vs.** Subparts



Examples

	Name ▲	Description	Modified	Organism	Sequenc...	Accession	Topology	Molecule Ty...	Genetic Code	Taxonomy	Path (Import...	Filename (Imported From)
	[1A]BBa_K314110	f1 origin	01 Jan 1980	.	455	BBa_K31411...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K314110-1A.gb
	[1B]BBa_K731722	T1 terminator from E. coli rrnB	01 Jan 1980	.	79	BBa_K73172...	linear	DNA	S			3.gb
	[1C]BBa_K398326	Promoter of the CaiF protein	01 Jan 1980	.	51	BBa_K39832...	linear	DNA	S			3.gb
	[1D]BBa_K731722	T1 terminator from E. coli rrnB	01 Jan 1980	.	79	BBa_K73172...	linear	DNA	S			3.gb
✓	[1E]BBa_K398331	pCaiF-B0032 measurement device	01 Jan 1980	.	935	BBa_K39833...	linear	DNA	S			3.gb
	[1F]BBa_K808025	FsC: Cutinase PET cleaving enzyme	01 Jan 1980	.	651	BBa_K80802...	linear	DNA	S			5.gb
	[1G]BBa_K314100	High constitutive expression cassette	01 Jan 1980	.	518	BBa_K31410...	linear	DNA	S			5.gb
	[1H]BBa_K808001	tphC: terephthalate periplasmatic binding prote...	01 Jan 1980	.	969	BBa_K80800...	linear	DNA	S			4.gb
	[1I]BBa_K314101	Low constitutive expression cassette	01 Jan 1980	.	518	BBa_K31410...	linear	DNA	S			1.gb
	[1J]BBa_K808003	tctB_162: small subunit B1 of the tripartite tric...	01 Jan 1980	.	489	BBa_K80800...	linear	DNA	S			1.gb
	[1K]BBa_K314201	Antitoxin Tsi2	01 Jan 1980	.	234	BBa_K31420...	linear	DNA	S			6.gb
	[1L]BBa_K808010	tphB: reaction from DCD to protocatechuate	01 Jan 1980	.	948	BBa_K80801...	linear	DNA	S			1.gb
	[1M]BBa_K314202	Toxin+ Antitoxin-Tse2 + Tsi2	01 Jan 1980	.	720	BBa_K31420...	linear	DNA	S			1.gb
	[1N]BBa_K808011	tphA1: reduces the tphA2A3 complex	01 Jan 1980	.	1,011	BBa_K80801...	linear	DNA	S			1.gb
	[1O]BBa_K346002	PmerT promoter (mercury-responsive)	01 Jan 1980	.	57	BBa_K34600...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K346002-1O.gb
	[1P]BBa_K808013	tphA3: Catalyzes together with tphA2 TPA to D...	01 Jan 1980	.	465	BBa_K80801...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K808013-1P.gb
	[2A]BBa_K548000	Humanized Aequorin	01 Jan 1980	.	594	BBa_K54800...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K548000-2A.gb
	[2B]BBa_K325210	Red Firefly Luciferase and LRE <i>i>L. Cruci...	01 Jan 1980	.	2,623	BBa_K32521...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K325210-2B.gb
	[2C]BBa_K648028	Cro, Lamda Repressor which activates the lytic...	01 Jan 1980	.	201	BBa_K64802...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K648028-2C.gb
	[2D]BBa_K325219	Red Firefly Luciferase and LRE (under pBAD)<...>	01 Jan 1980	.	3,841	BBa_K32521...	linear	DNA	Standard	-	C:\Users\Ce...	6411-BBa_K325219-2D.gb

iGEM Distribution 2021 0

Plate 1

380

D.gb

Plate 2

374

5.gb

Plate 3

348

4.gb

Plate 4

361

1.gb

Plate 5

380

1.gb


Plate 6

256

1.gb

Sequence View Annotations Text View Lineage Info

← → ↺ Extract R.C. Translate Add/Edit Annotation Allow Editing Annotate & Predict Save



10 % 🔍 🔍

Live Annotate & Predict

☐ Annotate from... ?

Source: Reference Features

Similarity: 100 %

Find: ☐ Best match ?

☒ All matching annotations

Apply Advanced...


☐ Find ORFs ?






Minimum size: 400

Genetic code: Standard

Start codon: CTG ATG TTG

	Name ▲	Description	Modified	Organism	Sequenc...	Accession	Topology	Molecule Ty...	Genetic Code	Taxonomy	Path (Import...	Filename (Imported From)
<input checked="" type="checkbox"/>	pSB1C3	High copy BioBrick assembly plasmid	29 Jul 2021 7:22 pm .		2,070	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	371808.gbk
<input checked="" type="checkbox"/>	pSB1C6	pSB1C3 based Type IIS backbone	29 Jul 2021 10:09 pm .		2,098	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	222895.gbk
<input checked="" type="checkbox"/>	pSB1K01	pOdd1 Loop Vector based on pSB1K3	30 Jul 2021 1:32 pm .		2,195	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	180932.gbk
<input checked="" type="checkbox"/>	pSB1K02	pOdd2 Loop Vector based on pSB1K3	30 Jul 2021 1:33 pm .		2,195	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	650996.gbk
<input checked="" type="checkbox"/>	pSB1K3	High copy BioBrick assembly plasmid	29 Jul 2021 10:18 pm .		2,204	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	174305.gbk
<input checked="" type="checkbox"/>	pSB1K03	pOdd3 Loop Vector based on pSB1K3	30 Jul 2021 1:34 pm .		2,195	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	761097.gbk
<input checked="" type="checkbox"/>	pSB1K04	pOdd4 Loop Vector based on pSB1K3	30 Jul 2021 1:35 pm .		2,195	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	695080.gbk
<input checked="" type="checkbox"/>	pSB1T3	High copy BioBrick assembly plasmid	29 Jul 2021 10:20 pm .		2,461	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	396566.gbk
<input checked="" type="checkbox"/>	pSB2K3	Inducible copy number BioBrick plasmid	29 Jul 2021 10:19 pm .		4,425	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	141419.gbk
<input checked="" type="checkbox"/>	pSB3C01	pEven1 Loop Vector based on pSB3C5	30 Jul 2021 1:28 pm .		2,729	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	735204.gbk
<input checked="" type="checkbox"/>	pSB3C02	pEven2 Loop Vector based on pSB3C5	30 Jul 2021 1:29 pm .		2,729	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	725183.gbk
<input checked="" type="checkbox"/>	pSB3C03	pEven3 Loop Vector based on pSB3C5	30 Jul 2021 1:30 pm .		2,729	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	127014.gbk
<input checked="" type="checkbox"/>	pSB3C04	pEven4 Loop Vector based on pSB3C5	30 Jul 2021 1:30 pm .		2,729	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	842518.gbk
<input checked="" type="checkbox"/>	pSB3C5	Low to medium copy BioBrick standard vector	29 Jul 2021 10:17 pm .		2,738	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	285332.gbk
<input checked="" type="checkbox"/>	pSB3K3	Low to medium copy BioBrick standard vector	29 Jul 2021 7:26 pm .		2,750	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	295979 (1).gbk
<input checked="" type="checkbox"/>	pSB3T5	Low to medium copy BioBrick standard vector	29 Jul 2021 10:21 pm .		3,252	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	814891.gbk
<input checked="" type="checkbox"/>	pSB4A5	Low copy BioBrick standard vector	29 Jul 2021 10:15 pm .		3,395	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	497914.gbk
<input checked="" type="checkbox"/>	pSB4C5	Low copy BioBrick standard vector	29 Jul 2021 10:18 pm .		3,221	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	475911.gbk
<input checked="" type="checkbox"/>	pSB4K5	Low copy BioBrick standard vector	29 Jul 2021 10:20 pm .		3,419	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	151165.gbk
<input checked="" type="checkbox"/>	pSB6A1	pBR322 converted to BioBrick vector (formally...	29 Jul 2021 10:16 pm .		4,022	<unknown	circular	DNA	Standard	-	C:\Users\Ce...	782671.gbk



 Hide fields	 Filter	 Group	 Sort								
<input type="checkbox"/>	Plate ▾	Well ▾	Part ▾	Type ▾	Short Description ▾	Plasmid ▾	Resistance ▾	Compatability ▾	Recognition Sites ▾	Sequencing ▾	Well Status ▾
41	Plate 1	3I	BBa_K538004	Coding	CspC (P. irgensii)	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
42	Plate 1	3J	BBa_S05060	Intermediate	K873000:B0015	pSB1C3	chloramphenicol	RFC[10] RFC[1000]	NotI	Confirmed(A)	OK
43	Plate 1	3K	BBa_K538000	Coding	Cpn10 (O. antarctica)	pSB1C3	chloramphenicol	RFC[10]	BsaI	Confirmed(A)	OK
44	Plate 1	3L	BBa_K774001	Regulatory	Mammalian-Bacterial Promoter: E9-ns2 CArG ...	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
45	Plate 1	3M	BBa_K525998	Regulatory	Promoter T7 and RBS	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
46	Plate 1	3N	BBa_K748002	Coding	Truncated lysostaphin coding sequence. ...	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
47	Plate 1	3O	BBa_K608002	Intermediate	strong Promoter and strong RBS	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
48	Plate 1	3P	BBa_S05048	Intermediate	B0034:K748000	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
49	Plate 1	4A	BBa_K331022	Composite	RBS with N-terminal Oligo Arginine - EYFP Fusion	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Bad Sequencing(A)	OK
50	Plate 1	4B	BBa_K381001	Composite	Nitrate reporter: PyeaR - GFP composite	pSB1C3	chloramphenicol	RFC[10]	BsaI	Confirmed(A)	OK
51	Plate 1	4C	BBa_K331024	Composite	RBS with N-terminal Oligo Arginine - ECFP Fusion	pSB1C3	chloramphenicol	RFC[10] RFC[1000]		Confirmed(A)	OK
52	Plate 1	4D	BBa_K314103	Composite	Lac induced expression cassette	pSB1C3	chloramphenicol	RFC[10] RFC[1000]	BsmBI BpiI	Long part(A)	OK
53	Plate 1	4E	BBa_K558009	Coding	Thiamin monophosphate pyrophosphorylase	pSB1C3	chloramphenicol	RFC[10]	BsaI	Confirmed(A)	OK
54	Plate 1	4F	BBa_K346000	Translational_Unit	RBS(B0032)+T3 DNA-directed RNA polymerase	pSB1C3	chloramphenicol	RFC[10]	BsaI BsmBI BpiI	Long part(A)	OK
55	Plate 1	4G	BBa_K575008	Composite	LasR/PAI1 Inducible	pSB1C3	chloramphenicol	RFC[10]	BsaI	Inconsistent(A)	OK
2,126 records											

Discussion